

## Simple protocol for Dig-RNA probes for in situ hybridization

(Wim Damen & Martin Klingler 12/99)

- Prepare plasmid (JETSTAR columns)
- Linearize plasmid (use enzyme that generates blunt end or 5' overhang)
- Purify by phenol/chloroform extraction and EtOH precipitation
- *In vitro* transcription:
  - 1  $\mu$ l DNA (500 ng)
  - 5  $\mu$ l H<sub>2</sub>O
  - 1  $\mu$ l 10x reaction buffer
  - 1  $\mu$ l rNTP mix **DIG** labeling (or rNTP mix **fluorescein** labeling)
  - 1  $\mu$ l RNase inhibitor (20 units)
  - 1  $\mu$ l RNA polymerase (T3/T7/Sp6)
  - 10  $\mu$ l
- Incubate 2 h at 37°C
- Add 90  $\mu$ l H<sub>2</sub>O, mix
- Run 5  $\mu$ l on a mini-gel (0.8% gel, 60 V, 15 min), place 25 ng template DNA in neighboring lane. RNA band should be 10 times stronger than template DNA.
  
- Add 1  $\mu$ l tRNA (20  $\mu$ g/ $\mu$ l)
- Add 16  $\mu$ l 3.2M LiCl, mix, and add 400  $\mu$ l EtOH
- Precipitate 2 h at -20°C
- Centrifuge 20 min at 15,000 rpm; pellet should be well visible
- Wash pellet with 500  $\mu$ l 70% EtOH (centrifuge again if pellet gets loose!)
- Dry in speedvac (about 5 min)
- Dissolve in 100  $\mu$ l H<sub>2</sub>O (DEPC-treated)
  
- Store aliquots at -80°C
- Test optimal probe concentration (try 1, 2, 4  $\mu$ l)